# MCScap

EEG accessories

---

**User Manual**

---

**EN**

EEG accessories MCScap. User manual.

Code: MCS.MC000000-63-EU-EN

Document revision: 5.2 of 2021-01-26

*\* The main version of the document is the English version.*

This document is a Manual to the use and maintenance of the EEG accessory MCScap and applies to the following accessories:

- electrode caps BASE, CLINIC, SLEEP, PROFESSIONAL;
- sets of textile caps and removable electrodes ELECTRODE SET;
- textile caps MCScap;
- electrodes MCScap-E / T / NT / NTC / NTH / C / TC.

*This document and Data Sheets\* for all accessories MCScap are available on [www.mcscap.com](http://www.mcscap.com) or upon request.*

*\* The Data Sheet for a specific accessory has a higher priority.*

## PRODUCT OFFICIAL INFORMATION

**Product name:** EEG accessories MCScap

**Description:** Accessories for non-invasive registration EEG and use of electroencephalographs and other amplifiers biological signals

**GMDN code:** 63131, Electroencephalographic electrode cap, reusable.

**Product code by 93/42/EEC:** MD 0100 General non-active, non-implantable medical devices

**Classification:** class I as per Appendix IX of Directive 93/42/EEC, rule 1.

**Declared for compliance with:** European Directives 93/42/EEC (MDD), 2011/65/EU (RoHS 2), 2002/96/EC (WEEE).

## CONTENT

---

<b>Symbols on labels</b> .....	<b>4</b>
<b>Warnings</b> .....	<b>5</b>
<b>Description</b> .....	<b>6</b>
Textile cap .....	6
Electrodes .....	8
Electrode caps .....	10
<b>Use</b> .....	<b>11</b>
<b>Cleaning</b> .....	<b>13</b>
<b>Disinfection and sterilization</b> .....	<b>14</b>
<b>Storage</b> .....	<b>14</b>
<b>Repairs</b> .....	<b>15</b>
<b>Utilization</b> .....	<b>15</b>
<b>Warranty</b> .....	<b>15</b>
<b>Appendix 1. Marking of the textile cap</b> .....	<b>16</b>
<b>Appendix 2. Main EEG electrode arrays</b> .....	<b>17</b>
<b>Appendix 3. Main pin layouts of common connector</b> .....	<b>19</b>

## SYMBOLS ON LABELS

 mks.ru	Consult instructions for use	 LOT	Batch code	 SN	Serial number	 REF	Catalogue number (article)
	Manufacturer	 EC REP	Authorized representative in the European Community		Separate collection for electrical and electronic equipment		Non sterile
	Not made with natural rubber latex		Keep away from sunlight		Temperature limit		Humidity limitation
	Wash by hand		Washing, mild process, maximum 30 Celsius		Do not tumble dry		Do not wring

## WARNINGS

- ⚠ Use the product only with compatible equipment. When choosing EEG accessories follow instructions for equipment with which they will be used.*
- ⚠ Use the product only for its intended purpose. Follow the recommendations of this document. When improper handling warranty may be canceled.*
- ⚠ It is prohibited to use the product in the working area of the magnetic resonance tomograph.*
- ⚠ Avoid using products near power installations and sources of high-frequency electromagnetic radiation (microwave ovens, x-ray machines, electrosurgical machines, mobile phones, etc.).*
- ⚠ Do not use product in water and in rain.*
- ⚠ Use products on different patients is allowed only after disinfection. Follow the instructions for cleaning and disinfecting products listed in this document.*
- ⚠ It is prohibited to impose the product on damaged areas of the skin.*
- ⚠ Do not overtighten the wires of the patient's body.*
- ⚠ Handle the components with care. Avoid excessive stretching of the material of the cap, twisting and tension of the wires.*
- ⚠ Avoid contact with the working surface of the electrodes with metals.*
- ⚠ Do not chlorinate electrodes.*

## DESCRIPTION

The MCScap accessories are designed for non-invasive registration EEG and use of electroencephalographs and other amplifiers biological signals. EEG accessories MCScap can be used in clinical neurology, functional diagnostics, and for scientific research.

### TEXTILE CAP

Textile cap is made of elastic material, preserving the shape and size. The cap provides the exact position of the electrodes on the head without additional measurements and adjustments. Large holes are provided for ventilation and access to the electrodes and patient's skin. The cap is fixed on the head with the chin or chest belt.

Textile cap with rings designed for fixing electrodes MCScap-E

Other models electrodes are installed directly into holes textile cap

Textile caps MCScap have identification size by color material or seam, and marking according to international system 10-20, 10-10 or 10-5\*.

\* Oostenveld, R. & Praamstra, P. *The five percent electrode system for high-resolution EEG and ERP measurements. Clinical Neurophysiology* 2001; 112: 713-719.

Table 1. MCScap size

Size	Head circumference	Color of cap for colored caps	Color of seams for black caps	Intended* for	Size of chest or chin belt
XL	60-66 cm	green 	 green	adults	adults
XL/L	57-63 cm	green and blue 	 green and blue	adults	
L	54-60 cm	blue 	 blue	adults(most)	
L/M	51-57 cm	blue and red 	 blue and red	teenagers, adults	children
M	48-54 cm	red 	 red	children, teenagers	
M/S	45-51 cm	red and yellow 	 red and yellow	children up to 5 years	
S	42-48 cm	yellow 	 yellow	children up to 2 years	
S/XS	39-45 cm	yellow and green 	 yellow and green	children up to 1 year	
XS	36-42 cm	green 	 green	babies up to 4 month	infant
Inf I	32-36 cm	blue 	 blue	infants up to 1 month	
Inf II	28-32 cm	red 	 red	newborn infants, preemies	
Inf III	24-28 cm	yellow 	 yellow	newborn infants, preemies	


\* Recommendations are based on WHO anthropological data <http://www.who.int/childgrowth> and ISO/TR 7250-2:2010. This information is a reference. In each case, the size of the cap should be selected, guided by the measured value of the circumference of the patient's head.

## ELECTRODES


Ag / AgCl electrodes MCScap guarantee minimal polarization and long-term stability of the signal.

*Ag/AgCl sintered* electrodes are produced according to technology solid sintering and are intended for special medical and scientific research.


*Ag/AgCl* electrodes are produced by technology surface deposition and are intended for routine clinical research.


MCScap-E


MCScap-T

MCScap-C  
MCScap-NTC

MCScap-TC  
MCScap-NT

MCScap-NTH

The conductive surface of electrodes does not have direct contact with skin, contact is provided with electrically conductive substance. A hole in the electrodes is provided to add a conductive gel.

For fastening electrodes on earlobes apply ear clamps.

Table 2. MCScap electrodes technical data

Model of electrode	MCScap-E	MCScap-T	MCScap-NT	MCScap-NTC	MCScap-NTH	MCScap-C	MCScap-TC
Material of conductive surface	Ag/AgCl sintered	Ag/AgCl sintered	Ag/AgCl sintered	Ag/AgCl sintered	Ag/AgCl sintered	Ag/AgCl	Ag/AgCl
Square of conductive surface	4 mm <sup>2</sup>	7 mm <sup>2</sup>	26 mm <sup>2</sup>	26 mm <sup>2</sup>	26 mm <sup>2</sup>	26 mm <sup>2</sup>	26 mm <sup>2</sup>
Internal diameter	8 mm	6 mm	8 mm	10 mm	7 mm	10 mm	8 mm
The square of contact of the electrode substance with the skin	50 mm <sup>2</sup>	28 mm <sup>2</sup>	50 mm <sup>2</sup>	79 mm <sup>2</sup>	39 mm <sup>2</sup>	79 mm <sup>2</sup>	50 mm <sup>2</sup>
Outer diameter	16 mm	12 mm	12 mm	15 mm	11 mm	15 mm	12 mm
The distance from the skin to the conductive surface	0.7 mm	1.6 mm	1.1 mm	3.2 mm	3.7 mm	3.5 mm	1.1 mm
The diameter of the hole in the electrode to add gel	2.8 mm						
Polarization	≤ 30 mV						
Resistance of electrodes insulation	≥ 1000 MΩ						
Dielectric strength of electrodes insulation	1500 V						
The impedance of the electrode	≤ 10 Ω						

## ELECTRODE CAPS

Electrode cap MCScap - this is textile cap with pre-installed Ag/AgCl electrodes and common connector or textile cap with set Ag/AgCl electrodes with individual connectors.


Table 3. MCScap electrode caps technical data \*

Model of textile cap	BASE	CLINIC	SLEEP	PROFESSIONAL	ELECTRODE SET
Electrode	MCScap-E	MCScap-C	MCScap-T	MCScap-NTC	MCScap-E
Number of electrodes	21				26
Connector type	Common connector DB-25M				Touch Proof 1.5 mm DIN 42 802-ST
Cable length	1.5 m				1.2 m
Mean time to failures (cycles **)	300	150	300	300	300

\* For custom products, these characteristics may be changed upon request.

\*\* Cycle means examining, washing and disinfecting the cap, cleaning and disinfecting the electrodes.


On the housing of the common connector contains the following information about the product:

The main schemes of electrode positions and common connector pin layouts are given in Appendices 2, 3.

## USE

1. Measure the circumference of the head of the patient and select the appropriate size for the cap.
2. Gently put on a cap on a patient's head starting from forehead and gradually pulling on his side nape.  
Electrode **Cz** should be located exactly in the middle of the line connecting bridge of the nose (*Nasion*) and occipital tubercle (*Inion*). The frontal electrodes **Fp1** and **Fp2** must be located directly above the superciliary arcs. Make sure that all the electrodes are located symmetrically about the middle line. If the location is incorrect, select a cap of a different size.
3. Fix the cap with a chin or chest belt.  
The tension of the belt must ensure reliable fixation of the cap. Avoid excessive strain. The patient should feel comfortable in the course of the entire study.  
 *When loose fit electrodes to head patient use cap smaller size or additional retainer.*

4. Plug in the cap connector to the amplifier and turn on the impedance control mode.

5. Treat your skin to reduce impedance.

In places where there are contact electrodes with skin, separate hair, and treat skin with special abrasive gel, following the instructions of the manufacturer.

 *Fat, lacquer for hair and other contamination skin increase impedance and impede quality registration EEG.*

6. Add an electrically conductive gel to each electrode.

The gel is injected with a syringe using a blunt needle through a hole in the electrode. In the process of adding the gel with the rotational movements of the needles, spread the hair under the electrode. Gel must be sufficient to provide electrical contact, but it does not have to go beyond the limits of the body of the electrode. Avoid connection stains electrically conductive gel between adjacent electrodes.

 *Use only special needles with blunt tip.*

 *It is recommended to use disposable needles. In the case of using reusable needles, before each use it is necessary to clean them and disinfect.*

## CLEANING

Electrodes and cap after use should be immediately cleared of gel while it has not dried up! A dried gel makes it difficult to achieve optimal electrode impedance and prevents quality recording of the EEG.

Electrode caps CLINIC, SLEEP, PROFESSIONAL recommended to be cleaned, not removing electrodes from textile cap. For cleaning electrode caps BASE and ELECTRODE SET it is recommended to remove electrodes MCScap-E from fixing rings.

 *Do not use for cleaning aggressive substances and shampoos.*

 *Avoid to put the connector in fluid!*

1. Disconnect the product from the amplifier.
2. Wash the product in water, using gentle cleaning agent and brush for cleaning electrodes.
3. Rinse the product under clean running water or distilled water.
4. Rinse with clean water and dry on a towel or hang.

*See video Cleaning of EEG accessories on [www.mcscap.com](http://www.mcscap.com).*

## DISINFECTION AND STERILIZATION

Disinfection of and sterilization of should be carried out only after procedure cleaning of the product.

For disinfection products it is recommended to use disinfectant solution, for example, Gigasept FF NEU or Bomix Plus. When choosing a disinfectant means consider *List of disinfectants VAH ([www.vah-online.de](http://www.vah-online.de))*. Please follow the instructions provided by the manufacturer of the selected agent. Acceptable use 70% solution ethanol method wipe.

-  *Do not soak the product in a disinfectant solution.*
-  *Do not use disinfectants with means of mass content shares of isopropyl alcohol more than 5%.*

Gas sterilization with ethylene oxide is acceptable.

-  *Do not use autoclave or other hot sterilization methods.*

## STORAGE

Temperature: from +5°C to +40°C. Humidity: from 20% to 80%.

-  *Don't expose direct sunlight.*

## REPAIRS

Textile cap and electrode user can replace his own.

For replacing defective electrode it is recommended to use MCScap repair kit.

*Instructions and video about replacing textile cap and electrode are available on [www.mcscap.com](http://www.mcscap.com)*

## UTILIZATION

This device must be disposed according to regulations of local legislation.


## WARRANTY

This warranty does not apply to any damages caused by wrong use, modification, and nonfulfillment of operation requirements.


The warranty period of storage in the manufacturer's packaging is **18 months** from the date of manufacture.

The warranty period for the product and reusable accessories is **six months** from the date of commissioning, but **not more than 24 months** from the date of manufacture. The warranty does not apply to consumable components.

## APPENDIX 1. MARKING OF THE TEXTILE CAP


Size XS, Inf I, Inf II, Inf III


Size XL-XS


## APPENDIX 2. MAIN EEG ELECTRODE ARRAYS


#39 (15 electrodes)

#23 (21 electrodes)


#40 (27 electrodes)


#41 (33 electrodes)


### APPENDIX 3. MAIN PIN LAYOUTS OF COMMON CONNECTOR


**#39-01 (15 electrodes)**

**#23-01 (21 electrodes)**

**#40-01 (27 electrodes)**

**#41-01 (33 electrodes)**

For other schemes, contact manufacturer.

**Manufacturer**

Medical Computer Systems Ltd.

Address: Passage 4922, bldg.4-2, Zelenograd, Moscow, 124460, Russia

Tel./fax: +7 (495) 913-31-94 / +7 (495) 913-31-95

E-mail: [mks@mks.ru](mailto:mks@mks.ru)

Internet: [www.mks.ru](http://www.mks.ru)

**European Authorized Representative**

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GVB-geliMED GmbH