MCScap

How to replace the electrode in
the cap

- CLINIC
- SLEEP
- PROFESSIONAL

EN

Instructions for replacing the electrode in the MCScap

Code: MCS.MC000000-67EN Revision: 1.3 of 29.03.2019

Purpose:

This instruction manual is intended for replacing the faulty electrode in MCScap of the CLINIC, SLEEP, PROFESSIONAL model.

Required tools and materials:

- Repair kit (CLINIC / SLEEP / PROFESSIONAL according to the MCScap model);
- Resistance measuring device (multimeter);
- Phillips screwdriver PH2;
- Diagonal pliers;
- A clerical or stripping knife (if necessary).

Before application this procedure, you have to sure that all requirement actions were done according to MCS. MC000000-63R (manual)

MCScap repair kit

- Repair electrode MCScap – C/T/NTC (electrode type according to MCScap model) – 2 pcs.
- Extractor – 1 pc.
- Cable tie 3*100 mm – 4/4/6 pcs.
- Cable fixator – 5 pcs.
- Thread 2 m – 1 pc.
- Instruction – 1 pc.

STEP 1 of 14

Identify the faulty electrode on the MCScap. Ensure that the electrode is defective with a multimeter. The device will show a contact break.

STEP 2 of 14

For silicone screeds and detachable fixators:
remove the cable from the defective electrode from the retainer. Be careful not to damage other electrodes!

For solid silicone tubes:
skip this step.


STEP 3 of 14

Remove the cable tie from the sleeve, bite the cable tie with pliers.


STEP 4 of 14

Disassemble the case of DB25 connector with a PH2 screwdriver.

STEP 5 of 14

Remove the cable tie from the DB25 connector, bite the cable tie with pliers.


STEP 6 of 14

Remove the crimped end of the cable from the connector using an extractor.

STEP 7 of 14


For a self-closing sleeves:

remove the cable of the defective electrode from the sleeve.

For non-self-closing sleeves:


1. Before removing, attach a 2 m thread to the end of the wire and cut the crimp;
2. Pull the electrode out of the braid from the side of the cap, while the attached thread must be inside the braid. The ends of the thread should be visible from both sides of the braid;
3. Unbind the thread from the cable of the defective electrode.


STEP 7.1 of 14

For solid silicone tubes:

remove the electrode out of the tubes and, if necessary, remove the tubes by cutting pliers.

STEP 8 of 14

For a self-closing sleeves:

Route the cable of the new electrode into the sleeve.


For non-self-closing sleeves:

1. Tie the crimped end of the wire of the new electrode to the thread from the side of the cap.


2. Extend the new electrode with a thread that is in the braid.

STEP 9 of 14

Insert the crimped end of the cable into the connector using an extractor.


STEP 10 of 14

Ensure that the installed electrode is working properly using a multimeter.

STEP 11 of 14

Install the cable tie in the DB25 connector.

STEP 12 of 14

Assemble the DB25 connector using a PH2 Phillips screwdriver.

STEP 13 of 14

Return the sleeve to the original place and install the cable tie on the sleeve.
Do not forget to pass the cable tie through the hole in the square retainer with velcro.

STEP 14 of 14

If necessary, install a retainer or cable tie to the cable of the new electrode, connecting it with the adjacent cables.

Fix the electrode in the hole of the cap in the place of the faulty one.

Your MCScap is ready!

Use the cap carefully, and it will last you a long time.

Thank you for choosing MCScap!

CONTACT INFORMATION

Medical Computer Systems Ltd.

Passage 4922, bldg.4-2, Zelenograd, Moscow, 124460, Russia

Tel.:+7 495 913-31-94.

E-mail: mks@mks.ru

Internet: www.mks.ru

Internet-shop MCScap: www.mcscap.com

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