

DATA SHEET

Adapter DB25-NeoRec16

Adapter for connecting the cap with a DB-25M common connector to an amplifier NeoRec 16.

REF

003-5-222

INTENDED USE

For connecting the electrode cap the amplifier.

SET

- Adapter DB25-NeoRec16.



DESCRIPTION

Adapter DB25-NeoRec16 is a passive adapter for connecting an EEG cap with a DB-25M common connector to the EEG amplifier NeoRec 16.

SPECIFICATION

Connector 1 type (to cap)	DB-25F
Connector 2 type (to amplifier)	ST40X-18S
Maximum number of electrodes to connect	17
Electrode positions	FP1, FP2, O1, O2, F3, F4, C3, C4, P3, P4, T7, T8, Fz, Cz, Pz, GND, AUX
Adapter / Extension lenght	0.3 ±0.01 m