

DATA SHEET

Chest belt MCScap

Belt for MCScap® fixing.

REF	Size
003-0-266	Adult
003-0-267	Children
003-0-268	Infant



INTENDED USE

For fixing EEG cap on the head.

SET

- Chest belt MCScap,
- User Manual.

DESCRIPTION

Elastic chest belt for fixing MCScap® caps with velcro for attaching to the cap.

Size of chest belt	Chest girth	Recommended for use with MCScap® sizes
Adult	70-120 cm	XL, XL/L, L, L/M
Children	40-70 cm	M, M/S, S, S/XS
Infant	30-40 cm	XS, InfI, InfII, InfIII

SPECIFICATION

Cap mount	velcro
Adjustment	for sizes "Adult", "Children"
Net weight	< 90 g (REF 003-0-266), < 55 g (REF 003-0-267), < 15 g (REF 003-0-268)
Gross weight	< 90 g (REF 003-0-266), < 55 g (REF 003-0-267), < 15 g (REF 003-0-268)