

EEG starter kit

Kit to start working with the EEG cap.

REF

003-7-001



INTENDED USE

To prepare for an EEG examination.

SET

- conductive gel,
- syringe,
- reusable needle,
- brushes for cleaning electrodes - 2 pcs.,
- spray for disinfection,
- operational documentation.

DESCRIPTION

EEG starter kit includes everything you need to start working with an EEG cap: contact electrode gel with high electrical conductivity, syringe, disposable needles for filling electrodes with gel through a special hole in the case of electrode, brushes for cleaning electrodes and spray for quick disinfection of electrodes and caps and accessories.

SPECIFICATION

Gel conductivity	3,0 - 4,0 S/m
Gel viscosity	9.0-11.0 Pa × s (shear rate 16.8±0.3 s ⁻¹); 17.0-21.0 Pa × s (shear rate 7.5±0.05 s ⁻¹)
Gel pH	6.0 - 7.0
Nominal gel weight	250 g
Nominal volume of syringe	25 ml
Needle length	15 mm
Needle diameter	2 mm

Nominal volume of the disinfectant	75 g
Net weight	405 ±5 g
Gross weight	410 ±5 g

