

DATA SHEET

Extension lead NeoRec16, 1.5 m

Extension lead for EEG amplifier NeoRec 16.



REF

003-5-223

INTENDED USE

For connecting the electrode cap the amplifier.

SET

• Extension lead NeoRec16, 1.5 m.

DESCRIPTION

The Extension lead NeoRec16, 1.5 m is designed to connect the EEG cap NeoRec cap to the EEG amplifier NeoRec 16 in case the amplifier is to be placed not on the cap, but in an arbitrary place on the human body or next to it.

SPECIFICATION

Connector 1 type (to cap)	ST60-18P
Connector 2 type (to ampifier)	ST40X-18S
Maximum number of electrodes to connect	17
Adapter / Extension lenght	1.5 ±0.1 m
Net weight	45 ±5 g