

## Set of adapters Snap-Banana for MCS ECG electrodes


Set of adapters for connecting MCS ECG electrodes to cables with a banana connector, 4 pcs. per pack.

### REF

007-0-305

### INTENDED USE

Connecting ECG electrodes to the amplifier.



### SET

- Snap-to-Banana adapter for MCS ECG electrodes – 4 pcs.

### DESCRIPTION

Snap-Banana adapters are passive adapters for connecting MCS ECG electrodes with a snap connector to lead cables with a banana connector.

### SPECIFICATION

Parameter	Value
Number of adapters	4 pcs.
Lead markings	presence
Type of connectors for electrodes	clamp for connector type snap
Connector for connecting to an amplifier	banana plug socket
Cable length	0.15 m

