

DATA SHEET

Adapter NeoRec21mini-DB25F

Adapter for connecting the cap with a DB-25M common connector to an amplifier NeoRec 21 mini.

REF

003-5-091

INTENDED USE

For connecting the electrode cap the amplifier.

SET

- Adapter NeoRec21mini-DB25F.



DESCRIPTION

Adapter NeoRec21mini-DB25F is a passive adapter for connecting an EEG cap with a DB-25M common connector to the EEG amplifier NeoRec 21 mini.

SPECIFICATION

Connector 1 type (to cap)	DB-25F
Connector 2 type (to amplifier)	Molex 52885-0374
Maximum number of electrodes to connect	22
Dimensions (LxWxH)	50 x 50 x 8 mm
Adapter / Extension lenght	-

