

DATA SHEET

## Soft skirt

Soft skirt to hold the electrically conductive gel and reduce the discomfort of skin contact with the hard electrode surface.

### REF

003-0-320



### INTENDED USE

To hold the EEG gel.

### SET

- Soft skirt - 10 pcs.

### DESCRIPTION

The Soft skirt is made of soft polyurethane and is designed to reduce the discomfort of skin contact with the hard electrode surface. Due to the additional volume of the gel, the skirt allows for electrical contact with a loose fit of the EEG electrode to the scalp. Soft skirt is designed for use with flat EEG electrodes MCScap-TC/T/NT/StTC/StNT.

### SPECIFICATION

Material	polyurethane
Dimensions (LxWxH)	13 x 13 x 5 ±1 mm
Weight of one item	< 0,5 g
Gross weight	< 5 g

